535 Personal Oxygen Station

DeVilbiss® iFill™ Personal Oxygen Station

Model 535D

CAUTION-Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING-Read and understand this instruction guide before operating this equipment.

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Improving People's Lives

TABLE OF CONTENTS

IEC Symbols	4
Important Safeguards	4
General Dangers and Warnings	4
Handling Warnings	5
Cautions & Notes	5
Provider Checklist	5
Introduction	5
Important Parts	
DeVilbiss iFill Personal Oxygen Station	6
iFill Oxygen Cylinder and Regulator	6
Accessory/Replacement Parts	6
Set-Up	
Choosing A Location	7
Transporting The Personal Oxygen Station	7
DeVilbiss iFill Personal Oxygen Station Operation Checklist	7
iFill Oxygen Cylinder External Examination	7
Operation	
Connecting The iFill Cylinder To The DeVilbiss iFill Personal Oxygen Station	8
Filling the iFill Oxygen Cylinder	8
Indicator Lights Explanation	8
Removing The iFill Oxygen Cylinder	9
Setting The iFill Oxygen Cylinder Rotary Selector To The Prescribed Setting	9
Troubleshooting	
DeVilbiss iFill Personal Oxygen Station	9
iFill Oxygen Cylinder/Regulator	10
Typical Questions and Answers	10
Cleaning/Maintenance	
Filter Cleaning	10
Cabinet	10
Specifications	
DeVilbiss iFill Personal Oxygen Station	11
iFill Oxygen Cylinder/Regulator	11
iFill Oxygen Cylinder Typical Fill Times	11

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IEC SYMBOLS

- | | |
|--|--|
|  Attention, consult instruction guide |  Type BF applied part |
|  AC current |  Double Insulated |
|  Electric Shock Hazard. Do Not Open |  Danger-No Smoking |

IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

- DANGER** Urgent safety information for hazards that will cause serious injury or death.
- WARNING** Important safety information for hazards that might cause serious injury.
- CAUTION** Information for preventing damage to the product.
- NOTE** Information to which you should pay special attention.

Important safeguards are indicated throughout this guide; pay special attention to all safety information.

READ ALL INSTRUCTIONS BEFORE USING.

SAVE THESE INSTRUCTIONS.

GENERAL DANGERS & WARNINGS

In order to ensure the safe installation, assembly, and operation of the personal oxygen station, these instructions **MUST** be followed:

DANGER-NO SMOKING

DANGER

- **DANGER-Electric Shock Hazard. Do not disassemble.** The DeVilbiss iFill Personal Oxygen Station contains no user serviceable parts. If service is required, contact your Sunrise Medical provider or authorized service center.
- **DANGER-Fire Hazard. Do not lubricate. Do not allow grease or oil from your hands or other source to come into contact with the regulator or cylinder valve connection.** These solutions may be flammable and cause injury.
- **DANGER-Keep the personal oxygen station at least 5 feet (1.6 m) from hot, sparking objects or naked sources of flame. Position your unit at least 6 inches (16 cm) from walls, draperies, or any other object that might prevent the proper flow of air in and out of your oxygen station. The oxygen station should be located so as to avoid pollutants or fumes.**
- **DANGER-Oxygen causes rapid burning. Do not smoke while your personal oxygen station is operating, or when you are near a person utilizing oxygen therapy. Keep the oxygen station at least 5 feet (1.6 m) from hot, sparking objects or naked sources of flame.**
- **DANGER-Do not lay the cannula down while the cylinder is delivering oxygen. High concentrations of oxygen can cause rapid burning.**

WARNING

- **Do not install, assemble, or operate this equipment without first reading and understanding this instruction manual and the oxygen cylinder warning label. If you are unable to understand the warnings, cautions and instructions, contact your provider or technical personnel if applicable before attempting to install or use this equipment - otherwise, injury or damage may occur.**
- **The use of this device is limited to an oxygen patient. Cylinders **MUST** be used **ONLY** by an oxygen patient and are not to be distributed to any other individual for any purpose.**
- **Do not use parts, accessories, or adapters other than those authorized by Sunrise Medical.**
- **NEVER block the air openings of the product or place it on a soft surface, such as a bed or couch, where the air opening may be blocked. Keep the openings free from lint, hair and the like.**
- **If the oxygen station has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or if it has been submersed in water, do not use and call a qualified technician for examination and repair.**

- If any leakage of the oxygen cylinder is detected, do not attempt to use the cylinder. Turn the rotary selector on the regulator to “OFF.” If leakage persists, place the cylinder outdoors and notify your Home Health Care Provider or Service Representative of this condition.
- The personal oxygen station is equipped with a High Pressure Relief Valve to ensure the user’s safety. When activated, this safety feature will make an extremely loud noise. If this noise occurs, turn the unit off and contact your Home Health Care Provider or Service Representative.
- Before moving or repositioning the personal oxygen station, always disconnect the AC power cord. Failure to do so may result in damage to the unit or personal injury.
- Children should always be supervised around the personal oxygen station. Failure to do so may result in damage to the unit or personal injury.
- Improper use of the power cord and plugs can cause a burn, fire, or other electric shock hazards. Do not use the unit if the power cord is damaged.
- For your safety, the oxygen cylinder must be used according to the prescription determined by your physician.
- Under certain circumstances, oxygen therapy can be hazardous. Seeking medical advice before using oxygen is advisable.

⚠ WARNING - HANDLING WARNINGS

- Use extreme care when handling and filling an oxygen cylinder. Full oxygen cylinders are under pressure and can become a projectile if dropped or mishandled.
- Never transport the oxygen station with a cylinder connected to it. Injury or damage can occur.

CAUTIONS & NOTES

CAUTION—Federal law restricts this device to sale by or on order of a physician.

NOTE— Do not connect to an electrical outlet controlled by a wall switch. No other appliances should be plugged into the wall outlet.

NOTE— The plug on the DeVilbiss iFill Personal Oxygen Station has one blade wider than the other. To reduce the risk of electric shock, this plug is intended to fit in a wall outlet only one way. Do not attempt to defeat this safety feature.

NOTE— Sunrise Medical recommends leaving a full cylinder of oxygen with the patient **AFTER** setting up the personal oxygen station.

PROVIDER CHECKLIST

Before leaving a DeVilbiss iFill Personal Oxygen Station with a user, the following checklist **MUST** be completed:

1. Check all parts for shipping damage. In case of damage, do not use. Contact Sunrise Medical for further instructions.
2. Make sure the hydrostatic test on the oxygen cylinder has not expired.
3. Instruct the user on the safe operation of the personal oxygen station and review ALL warnings.
4. Leave a copy of this instruction guide with the user.

NOTE— Sunrise Medical recommends leaving a full cylinder of oxygen with the patient **AFTER** setting up the personal oxygen station, and instructing the patient to always keep a filled cylinder on hand.

INTRODUCTION

This instruction guide will acquaint you with the DeVilbiss iFill Personal Oxygen Station. Make sure that you read and understand this guide before operating your oxygen station. Important safeguards are indicated throughout this guide; pay special attention to all safety information. Contact your Sunrise Medical/DeVilbiss equipment provider should you have any questions. **NOTE**– All references to cylinder throughout the guide represent an oxygen cylinder.

INDICATIONS FOR USE

The Model 535D Oxygen Cylinder Filling System is intended for use in supplying pressurized oxygen to fill oxygen cylinders for patients' ambulatory use. The device is intended to provide 93% (±3%) oxygen. This device can be used in homes, nursing homes, patient care facilities, etc.

IMPORTANT PARTS

DEVILBISS IFILL PERSONAL OXYGEN STATION (FIGURES A & B)

1. Oxygen Fill Connector (referred to as Fill Connector throughout guide)
2. Oxygen Fill Connector Cover/Cylinder Adapter (referred to as Fill Connector Cover and/or Cylinder Adapter throughout guide)
3. Oxygen Fill Port Collar (referred to as Fill Port Collar throughout guide)
4. Instruction Label
5. Indicator Panel Label
 - a. Service Required (Red) LED
 - b. Filling (Green) LED
 - c. Full (Green) light
 - d. Standby (Green) light
6. Start/Stop Switch
7. Cradle
8. Circuit Breaker
9. Front Handle
10. Rear Handle
11. Carry Handle Recess
12. Intake Filter(s) (2 - one on each side)
13. Hour Meter - behind filter on Indicator Panel label side (right side)


Figure A

IFILL OXYGEN CYLINDER & REGULATOR (FIGURE C)

1. Rotary Selector
2. Nipple Connector
3. Hydrostatic Test Date - Contact your homecare provider or the cylinder manufacturer for details.
4. Cannula Fitting


Figure B

ACCESSORY/REPLACEMENT PARTS

1. Intake Filter(s)
2. iFill Oxygen Cylinder w/Integrated PulseDose Regulator (US)

- 535D-605 (6 filters)
- PD1000A-M4
- PD1000A-M6
- PD1000A-ML6
- PD1000A-C
- PD1000A-D
- PD1000A-E
- PD1000A-CAN-M4
- PD1000A-CAN-M6
- PD1000A-CAN-ML6
- PD1000A-CAN-C
- PD1000A-CAN-D
- PD1000A-CAN-E

iFill Oxygen Cylinder w/Integrated PulseDose Regulator (CAN)

3. iFill Oxygen Cylinder w/Continuous Flow Regulator (US)

- 535D-M6-CF
- 535D-ML6-CF
- 535D-C-CF
- 535D-D-CF
- 535D-E-CF
- 535C-M6-CF
- 535C-ML6-CF
- 535C-C-CF
- 535C-D-CF
- 535C-E-CF

iFill Oxygen Cylinder w/Continuous Flow Regulator (CAN)


Figure C

SET-UP

NOTE– Before operating the DeVilbiss iFill Personal Oxygen Station for the first time, locate the device where desired.

CHOOSING A LOCATION

⚠ WARNING

The personal oxygen station should be at least six (6) inches away from walls, draperies, or furniture to ensure sufficient airflow.

Avoid placing the oxygen station next to heaters, radiators, or hot air registers. It should be kept at least 5 feet (1.6m) away from hot, sparking objects or naked sources of flame.

Do not use in a closet.

Do not use in an area where the air may be contaminated with carbon monoxide as this may shorten the life of the unit (i.e. near running gasoline engines, furnace, or heater).

TRANSPORTING THE DEVILBISS IFILL PERSONAL OXYGEN STATION

⚠ WARNING

NEVER transport the Personal Oxygen Station with a cylinder connected to or stored on it, otherwise injury or damage can occur.

CAUTION– When transporting the DeVilbiss iFill Personal Oxygen Station, be careful not to bump the unit or the connectors into obstacles. Otherwise damage to the oxygen station may occur.

There are three (3) methods to safely transport the DeVilbiss iFill Personal Oxygen Station (Figure 1):

1. Pull or push the unit on its casters using the front handle. Travel in a line oriented front to back to avoid tipping. Stop the unit to change direction.
2. Carry the unit using the front and rear handles simultaneously.
3. Carry the unit using the front handle and the recess located under the rear of the base.


Figure 1

DEVILBISS IFILL PERSONAL OXYGEN STATION OPERATION CHECKLIST

NOTE– Each time the Oxygen Station is used to fill a cylinder, complete the following steps:

1. Check the Hydrostatic Test Date **Figure C, 3** - Aluminum cylinders must undergo testing every five (5) years. Contact your homecare provider or the cylinder manufacturer for details.

NOTE– Do not fill oxygen cylinders that have not been tested in the past five (5) years. Contact your Sunrise provider for replacement.

2. Perform the prefill inspection on the cylinder by following *iFill Oxygen Cylinder External Examination* below.

⚠ WARNING

ONLY use cylinders that have the DeVilbiss iFill Personal Oxygen Station fill connection.

All cylinders must be inspected before attempting to fill - otherwise, injury or damage may occur.

IFILL OXYGEN CYLINDER EXTERNAL EXAMINATION

1. Examine the outside of the cylinder for the following conditions, and replace the cylinder if they exist:
 - a) Dents or dings
 - b) Arc Burns
 - c) Oil or Grease
 - d) Any other signs of damage that might cause a cylinder to be unacceptable or unsafe for use.
2. Examine the cylinder for evidence of fire or thermal damage. Evidence includes charring or blistering of the paint, or other protective coating or heat sensitive indicator. If fire or thermal damage is found, replace the cylinder.

3. Inspect the oxygen fill connector for the following:
 - a) Debris, oil or grease
 - b) Noticeable signs of damage
 - c) Signs of corrosion inside the valve
 - d) Signs of excessive heat or fire damage
- Do not use oxygen cylinder if found, call your homecare provider.

OPERATION

CONNECTING THE IFILL OXYGEN CYLINDER TO THE DEVILBISS IFILL PERSONAL OXYGEN STATION

⚠ DANGER

Fire Hazard. Do not lubricate. Do not allow grease or oil from your hands or other source to come into contact with the regulator or cylinder valve connection. These solutions may be flammable and cause injury.

⚠ WARNING

Do not modify ANY connections on the personal oxygen station. NEVER use tools of any kind to connect/disconnect the cylinder and the oxygen station. Severe injury and/or damage may occur.

Do not drop oxygen cylinders. Use two (2) hands when handling/transporting oxygen cylinders. Otherwise, injury or damage may occur.

NOTE– Refer to the *GENERAL DANGERS & WARNINGS* and the *HANDLING WARNINGS* in this manual.

1. Plug the power cord into the wall outlet.
 2. Set the rotary selector on the cylinder to “OFF.”
 3. Remove the oxygen fill connector cover from the fill connector.
 4. Position the cylinder over the cradle while aligning the nipple connector on the oxygen cylinder with the fill connector. Press until cylinder “Clicks” into place (Figure 2).
- NOTE**– When using the smaller M4 or M6 size cylinders, it is necessary to place the oxygen fill connector cover/cylinder adapter in the cradle to help support the cylinder. When using the larger ML6, C, D, or E cylinders, the oxygen fill connector cover/cylinder adapter is not necessary and can be allowed to hang over the side.


Figure 2

FILLING THE IFILL OXYGEN CYLINDER

NOTE– Refer to the *GENERAL DANGERS & WARNINGS* and the *HANDLING WARNINGS* in this manual.

1. Make sure the personal oxygen station is plugged in.
2. Make sure the cylinder is securely connected to the oxygen station. Refer to *Connecting The iFill Oxygen Cylinder To The DeVilbiss iFill Personal Oxygen Station*.
3. Push the start/stop switch on the control panel (Figure 3).
4. The following sequence of events should occur:
 - a. The FILLING (green) light will illuminate while the cylinder is filling (refer to *iFill Oxygen Cylinder Fill Times*).
 - b. The FULL (green) light will illuminate when the cylinder is full. Proceed to *Removing the iFill Oxygen Cylinder*.

NOTE– If the Service (RED) light and audible alert are on, proceed to *TROUBLESHOOTING* in this instruction manual.


Figure 3

INDICATOR LIGHT EXPLANATION

INDICATOR LIGHT COLOR	OXYGEN STATION STATUS	EFFECT
None	Oxygen station is unplugged.	None
Standby (green)	Oxygen station is ready to begin filling cylinders.	Cylinder filling not started
Filling (green)	Oxygen station is on and cylinder is filling.	Cylinder is filling
Full (green)	Cylinder is full.	Remove cylinder.
Service (red) and Audible Alert	Oxygen station is on but cylinder is not filling due to internal failure.	Contact your Sunrise Medical Provider

REMOVING THE IFILL OXYGEN CYLINDER

⚠ WARNING

Do not modify any connections on the personal oxygen station. Never lubricate the connections. Never use tools of any kind to connect/disconnect the oxygen cylinder and the oxygen station. Otherwise, severe injury and/or damage may occur.

Use extreme care when handling and filling an oxygen cylinder. Full oxygen cylinders are under pressure and can become a projectile if dropped or mishandled.

NOTE– Refer to the *GENERAL DANGERS & WARNINGS* and the *HANDLING WARNINGS* in this manual.

1. Remove the full oxygen cylinder by pushing down on the fill port collar located under the cylinder nipple connector while using your other hand to steady the cylinder (Figure 4).
2. Lift up on the cylinder to remove from the oxygen fill connector. The green Standby light will turn on after a short delay.
3. **Oxygen will escape from the fill connector for a short time after the cylinder has been removed. This may be accompanied by a “pop” with a rush of air. This is normal.**
4. Perform one (1) of the following steps:
 - A. If desired, fill another oxygen cylinder. Refer to the *DeVilbiss iFill Personal Oxygen Station Operation Checklist*.
 - B. Place the oxygen fill connector cover onto the fill connector.


Figure 4

NOTE– The oxygen fill connector cover should always be used whenever cylinders are not being filled.

SETTING THE IFILL OXYGEN CYLINDER ROTARY SELECTOR TO THE PRESCRIBED SETTING

1. Attach the nasal cannula to the cannula fitting of the iFill oxygen cylinder.
2. Turn the rotary selector to the L/min. setting prescribed by your physician.

⚠ WARNING

Changing the L/min. setting on the rotary selector will affect the dose of oxygen delivered, DO NOT readjust the L/min. setting unless directed by your physician.

3. To turn the flow of oxygen off, turn the rotary selector counterclockwise to the “OFF” position.

TROUBLESHOOTING

The following troubleshooting chart will help you analyze and correct minor system malfunctions. If the suggested procedures do not help, call your Sunrise Medical/DeVilbiss homecare provider. Do not attempt any other maintenance.

⚠ WARNING

To avoid electric shock hazard, do not remove the cabinet. The cabinet should only be removed by a qualified Sunrise Medical/DeVilbiss homecare technician.

DEVILBISS IFILL PERSONAL OXYGEN STATION

PROBLEM	SOLUTION
No indicator lights on.	1. Check that the personal oxygen station power cord is properly plugged into the wall outlet. Ensure the circuit breaker near the power cord entry is pressed in. If there are still no indicator lights, contact your Sunrise Medical/DeVilbiss provider.
Red light on.	1. Ensure the iFill oxygen cylinder is connected correctly. 2. Press the start/stop switch. If the Red light remains on, contact your Sunrise Medical/DeVilbiss provider for service. 3. If the Red light turns on during the filling process, check that the cylinder is connected properly. Remove the cylinder and reinstall. Ensure the cylinder rotary selector is in the OFF position. Press the start/stop switch to turn off the Red light and then again to restart the cylinder filling process. If the Red light turns on again during the filling process, contact your Sunrise Medical/DeVilbiss provider for service.
Excessive fill times OR Green filling light stays on.	1. Ensure the iFill oxygen cylinder rotary selector is set to OFF. 2. Check the fill times shown in <i>Typical Oxygen Cylinder Fill Times</i> chart. 3. Check that the cylinder is connected properly: Press the start/stop switch to return to Standby. Remove the cylinder and reinstall. Press the start/stop switch to begin filling the cylinder. 4. If problem persists, contact your Sunrise Medical/DeVilbiss provider for service.

IFILL OXYGEN CYLINDER/REGULATOR

PROBLEM	SOLUTION
Hissing sound from iFill oxygen cylinder	<ol style="list-style-type: none"> 1. Turn rotary selector OFF. 2. If hissing sound continues, call your Sunrise Medical/DeVilbiss provider.
iFill oxygen cylinder is on but oxygen is not being delivered.	<ol style="list-style-type: none"> 1. Refill cylinder. 2. Check that cannula is straight and not pinched. If damaged, replace. 3. Call your Sunrise Medical/DeVilbiss provider for service.

TYPICAL QUESTIONS AND ANSWERS

- Q.** *Why does my iFill make a “pop” sound when I remove the cylinder?*
- A.** The “pop” sound is pressure being released from the system to make the next cylinder easier to install.
- Q.** *Why is the coupler only to be cleaned with a clean, dry, lint-free cloth and only when necessary?*
- A.** Cleaning solutions may be flammable and/or may contaminate the cylinder.
- Q.** *Why can't I use my old cylinders with the iFill?*
- A.** The iFill has a proprietary connection as required by the FDA and, therefore, other cylinders are not compatible.

CLEANING/MAINTENANCE**DEVILBISS IFILL PERSONAL OXYGEN STATION FILTER**

CAUTION – *Do not operate the personal oxygen station without the filters installed.*

NOTE– *The DeVilbiss iFill Personal Oxygen Station has two (2) filters, one on each side of the cabinet (FIGURE A).*

NOTE– *Perform this procedure at least ONCE A WEEK.*

1. Unplug the personal oxygen station from the wall outlet.
2. Remove the filters.
3. Clean the filters with a vacuum cleaner or wash in warm soapy water and rinse thoroughly.
4. Dry the filters thoroughly before reinstallation.

EXTERIOR CABINET** WARNING**

Electric Shock Hazard: Unplug the DeVilbiss iFill Personal Oxygen Station when cleaning. Do not remove the oxygen station cabinet. The cabinet should only be removed by a qualified Sunrise Medical/DeVilbiss technician. Do not apply liquid directly to the cabinet nor utilize any petroleum-based solvents or cleaning agents.

1. Clean the exterior cabinet by using a damp cloth or sponge with a mild household cleaner and wipe it dry.
2. Only if necessary, wipe the coupler with a clean, dry, lint-free cloth.

SPECIFICATIONS

DEVILBISS IFILL PERSONAL OXYGEN STATION

Operating Temperature:41° to 95°F (5° to 35°C)
 Operating Humidity Range:15 to 95% R.H. non-condensing
 Operating Altitude:0 to 6,562 Feet (0 to 2,000 Meters)
 Storage Temperature Range:-40°F to +158°F (-40°C to +70°C)
 Storage Humidity Range:10 to 95% non-condensing
 Electrical Rating:115V~ 60Hz 4.1 Amps
 Operating Voltage Range:97 – 127V~ 60HZ
 Power Range:400 Watts Average
 Oxygen Purity:93% O2 ± 3%
 Dimensions: (including casters) w/o cylinder12.25"W x 28.5"H x 22.5"D
 Weight:66 lbs
 Shipping Weight:76 lbs
 This unit is classified as nonprotected (ordinary equipment) per EN60601-1IPX0

Specifications subject to change without notice.

IFILL OXYGEN CYLINDER TYPICAL FILL TIMES

NOTE- All filling times are approximate and may vary depending on altitude and environmental conditions.

Typical cylinder fill time from empty to 2,000 +/- 200 psig are as follows:

- M4 60 Minutes
- M6 75 Minutes
- ML6 90 Minutes
- C 130 Minutes
- D 215 Minutes
- E 350 Minutes

NOTE- Degradation of performance may occur if unit is operated outside of specified operating parameters.

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Sunrise Medical Pty. Limited