# FLOW-SAFE II EZ®

DISPOSABLE CPAP SYSTEM

**DESCRIPTION:**  
The Mercury Flow-Safe II EZ Continuous Positive Airway Pressure (CPAP) device is a respiratory aid intended for use with a face mask, nebulizer and gas supplying device to elevate pressure in the patient's lungs while delivering aerosolized medication.

**CONNECTIONS:**

- Standard oxygen tubing nipple
- Patient connection ISO 5356-1 - 15mm taper female and 22mm taper male
- Nebulizer connection 18mm taper female and 22mm taper male

**INDICATIONS FOR USE:**

The Mercury Flow-Safe II EZ CPAP device is intended to provide CPAP to spontaneously breathing patients in the hospital and pre-hospital environment.

**WARNINGS:**

- Do not allow smoking or use near smoking equipment, open flame, oil or other flammable chemicals.

**CAUTIONS:**


- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Do not clean, soak, rinse or sterilize.
- Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.
- In the event of undesirable flow rate from oxygen source, simply remove the device and place on supplemental oxygen per protocol.
- Use of the Flow-Safe with non-back pressure compensated flow devices may affect input gas liter flow. Always verify delivered CPAP pressure on manometer.
- Activation or deactivation of nebulizer may affect the delivered CPAP pressure. Always verify CPAP pressure with manometer.
- Flowmeters capable of delivering up to 25 LPM may be required to operate both CPAP device and nebulizer simultaneously.
- Use of any nebulizer other than the one supplied may affect performance.

**CONTRAINdications:**

- Respiratory Arrest
- Unconscious
- Cardiogenic Shock
- Pneumothorax
- Facial Anomalies
- Facial Trauma

**DIRECTIONS FOR USE:**

1. Connect O<sub>2</sub> tubing nipple to gas source.
2. Place the face mask securely to the patient's face using head harness.
3. With nebulizer in the OFF position (**Figure 1**), slowly increase gas flow to 6 or 8 LPM. Check face mask fit to patient and device connection for leaks.
4. Adjust the device until desired pressure is obtained. Flow of 12 - 14 LPM is required to reach CPAP pressure of 8.5 - 10 cm H<sub>2</sub>O.
5. Do not exceed 30 LPM.
6. Patient SaO<sub>2</sub> should be monitored using a pulse oximeter.
7. To activate nebulizer rotate knob to ON position (**Figure 2**).
8. If necessary, readjust flowmeter to obtain desired CPAP pressure. Up to 25 LPM may be required.

**MEASURING PRESSURE:**  
- Pressure relief limits maximum CPAP pressure to 25 cm H<sub>2</sub>O @ 25 LPM.  
- Do not exceed pressure limit of manometer (25 cm H<sub>2</sub>O).  
- Manometer accuracy  $\pm 3$  cm H<sub>2</sub>O up to 15 cm H<sub>2</sub>O and  $\pm 5$  cm H<sub>2</sub>O over 15 cm H<sub>2</sub>O.

**SPECIFICATIONS:**

Sample Guidelines for Preparing RX dosing

FLOWMETER SETTING l/min	14 ~ 15	23 ~ 24
CAP AP PRESSURE cm H <sub>2</sub> O	4 ~ 5	9 ~ 10
FLOW THROUGH EZFLOW MAX	6 l/min	10 l/min
OUTPUT	12 mL/hour	16 mL/hour
RX (mg/hour)	5 10 15 20	5 10 15 20
Treatment Duration (hours)	1 2 1 2 1 2 1 2 1 1.5 1 1.5 1 1.5	1 2 1 2 1 2 1 2 1 1.5 1 1.5 1 1.5
Medication @ 5 mg/mL (mL)	1 2 2 4 3 6 4 8 1 1.5 2 3 3 4.5 4 6	1 2 2 4 3 6 4 8 1 1.5 2 3 3 4.5 4 6
Saline (mL)*	11 22 10 20 9 18 8 16 15 22 14 21 13 20 12 18	11 22 10 20 9 18 8 16 15 22 14 21 13 20 12 18

\* Rounded to the nearest mL

**EXAMPLE:**  
To deliver 10 mg/hour of 5 mg/mL medication at 10 L/min; 1.5 Hour Duration: Place 3 mL of medication + 21 mL of Saline into nebulizer reservoir and run for 1.5 hours.

**CAUTION:** The listed output is nominal value only, actual output may vary depending on device, accuracy of flowmeter, and the existence of tubing connection leaks. Follow up output checks are recommended and flowmeter adjustments may be necessary.

© Registered U.S. Trademark of Mercury Enterprises, Inc.  
Patent # US 8,522,618 B1, US 6,338,443 B1, US 9,370,635 B2, US 9,511,202, Other Patents Pending

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