

Sunflower Warmer

Radiant Warmer



Simple and Versatile

Total support from postnatal treatment to NICU care.
The most suitable model can be chosen from a cost-effective line up.

Hi-Low/Tilt



Cabinet/Tilt



Cabinet



Wall Mount



Stand



*Sufficient optional accessories to neonatal care
can be mounted.*




Hi-Low/Tilt

* The drawer shown is
optional

Surgical operations can be performed safely in the NICU.

◆ **Wide accessible workspace for surgical treatment**


- Easy to handle the anesthetic equipment.
- A range of peripheral devices such as resuscitation unit (T-piece resuscitator, oxygen blender and suction unit), phototherapy unit, examination light, weight monitor, syringe pump, etc. can be mounted through F-Rail.




- LED light helps the medical staff to observe the infant.

*“Baby-friendly” function
to relieve infant’s stress.*


◆ **Reduces the infant's stress
and supports daily care giving**



- Baby-guard with very quiet rotary damper.



- Double-wheel caster enables stable and smooth transportation.



- Self closing doors that close slowly and silently.

*The drawer shown is optional.




Best for urgent postnatal treatment.

Prompt and safe resuscitation can be performed.

◆ **Resuscitation treatment in compliance with ILCOR Consensus 2010 can be achieved.**

• CPR Timer prompts the user to assess the health condition of the infant every 30 seconds.



• Resuscitation unit (T-piece resuscitator, oxygen blender and suction unit) and weight monitor, etc. can be mounted to F-Rail.



Cabinet/Tilt



Utility design.

Sunflower Warmer
Hi-Low/Tilt



HL Stand	✓
Fixed Stand	—
Apgar Timer	✓
CPR Timer	✓
Mattress platform tilting	✓
Rotary damper	✓
Drawer	▲
Cabinet	—
Double-wheel Caster	✓
Single-wheel Caster	—
X-Ray cassette tray port	✓

✓ : Included as standard
— : Not available ▲ : Optional

Sunflower Warmer
Cabinet/Tilt



HL Stand	—
Fixed Stand	✓
Apgar Timer	✓
CPR Timer	✓
Mattress platform tilting	✓
Rotary damper	✓
Drawer	—
Cabinet	✓
Double-wheel Caster	✓
Single-wheel Caster	—
X-Ray cassette tray port	✓

✓ : Included as standard — : Not available

Sunflower Warmer
Cabinet



Steel



Wood

HL Stand	—
Fixed Stand	✓
Apgar Timer	✓
CPR Timer	✓
Mattress platform tilting	—
Rotary damper	✓
Drawer	—
Cabinet	✓
Double-wheel Caster	✓
Single-wheel Caster	—
X-Ray cassette tray port	✓

✓ : Included as standard — : Not available

Maneuverable design.

Sunflower Warmer
Stand



HL Stand	—
Fixed Stand	✓
Apgar Timer	✓
CPR Timer	✓
Mattress platform tilting	—
Rotary damper	—
Drawer	—
Cabinet	—
Double-wheel Caster	—
Single-wheel Caster	✓
X-Ray cassette tray port	—

✓ : Included as standard — : Not available

Space-saving design.

Sunflower Warmer
Wall Mount



HL Stand	—
Fixed Stand	—
Apgar Timer	✓
CPR Timer	✓
Mattress platform tilting	—
Rotary damper	—
Drawer	—
Cabinet	—
Double-wheel Caster	—
Single-wheel Caster	—
X-Ray cassette tray port	—

✓ : Included as standard — : Not available







Atom's code. 60983
Drawer, 2 Layers
For Infa Warmer i/
Sunflower Warmer
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
Self-closing drawer that
closes slowly and silently.



Atom's code. 60899
Handle
For F-Rail
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
The vertical position of the
handle can be adjusted
easily according to the
height of the operator.
It is conveniently used to
move Sunflower Warmer to
another place.



Atom's code. 61616
MF-Rail
For F-Rail
(Common to Left & Right)
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
This rail is used to mount the
peripheral devices attachable
to the MF-Rail for F-Rail.



Atom's code. 60893
I.V Pole
For F-Rail
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
With the load capacity of
10kg, this I.V Pole can mount
the weight of various
peripheral devices.



Atom's code. 61642
Drainage Hanger
For Infa Warmer i/
Sunflower Warmer
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
The Drainage Hanger is used
to hang a low-pressure
continuous aspirator and
drainage bag on Sunflower
Warmer.



Atom's code. 60871
Tube Holder
For Infa Warmer i/
Sunflower Warmer/V-505
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
The Tube Holder is attached
to the baby guard at the head
end of the warmer to hold
the patient circuit securely in
place.



Atom's code. 61647/61640
(for USA only)
Weight Monitor Unit
For Infa Warmer i/
Sunflower Warmer
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
You can weigh a baby
placed on the mattress
platform whenever necessary.



Atom's code. 60216
InfaLite
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
With a LED light and hands-free
motion sensor, both brightness
and hygiene have been
considered with this
examination light.

*Either F-Rail attachment or Pole
attachment is required for the
installation of InfaLite.



Atom's code. 60984
Shelf
For Sunflower Warmer
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
Shelf with wide space can
place a SpO2 monitor, etc.



Atom's code. 61618
Tray Set (S)
For F-Rail
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
The vertical position of the
Tray Set (S) can be adjusted
easily according to the height
of the operator.
Small instruments necessary
for routine medical
procedures can be placed
on this small tray.



Atom's code. 60362
Resusci Flow
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
"Resusci Flow" is a resuscitation unit that performs
artificial respiration with a stable pressure.

•Simple operation supports resuscitation.
•Safety function to protect newborn's respiration.
•Supports artificial respiration in various situations.



Atom's code. 60894
Basket
w/arm For F-Rail
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
The vertical position of the
basket can be adjusted
easily as desired.
The basket is deep enough
to conveniently accommodate
tall containers such as bottles
of alcohol disinfectant
solutions.



Atom's code. 60350
Monitor Table
for Philips monitor
For F-Rail
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
Philips monitor can be
mounted on the monitor table.



Atom's code. 60361
BILI-THERAPY Spot Type Arm Type
Hi-Low/Tilt Cabinet/Tilt Cabinet
Stand Wall Mount
Effective Phototherapy Treatment and Space-saving
•Effective phototherapy treatment by blue LEDs
Promote early hospital discharge by reducing the treatment time in providing high irradiation
intensity and effective wavelength.
•Space-saving
The unit can be mounted on the F-Rail without using a stand for space-saving
around the incubator or warmer.
Switchable to examination lamp with one-touch of a button
The condition of an infant can be monitored by the examination lamp (White LEDs).
The dimmer function ideal for developmental care is available.

Products Lineup

Sunflower Warmer		
Type	Atom's code	Description
Hi-Low/Tilt	61564	AC120V±10%
	61565	AC230V±10%
Cabinet/Tilt	61566	AC120V±10%
	61567	AC230V±10%
Cabinet/Steel	61568	AC120V±10%
	61569	AC230V±10%

Type	Atom's code	Description
Cabinet/Wood	61570	AC120V±10%
	61571	AC230V±10%
Stand	61572	AC120V±10%
	61573	AC230V±10%
Wall Mount	61574	AC120V±10%
	61575	AC230V±10%

Manufacturer: ATOM MEDICAL CORPORATION

3-18-15, Hongo, Bunkyo-ku, Tokyo, JAPAN
<http://www.atomed.co.jp>

Exporter: ATOM MEDICAL INTERNATIONAL, INC.

Iwakata Bldg., 3-18-16, Hongo, Bunkyo-ku, Tokyo, JAPAN
Phone: + 81 3 3815 2941 Fax: + 81 3 3812 9670/4080